

High Noble, Noble, Non-Precious *with*



Griffin Dental Laboratories produces full cast restorations because of their strength, durability and bio compatibility.

Prep Requirements

- Minimal Prep requirements. Can be used in areas with minimal occlusal reduction.

Indications

- These restorations can be used for a wide variety of cases as listed below
 - Crowns
 - Bridges
 - Veneers
 - Inlays
 - Onlays

Contraindications

- Alloys not to be used in the fabrication of
 - Partial(RPDs) or
 - Implants

Cementation

Conventional Cementation is used in all instances. Follow manufactures recommendations.

Alloy Types

High Noble, Noble, Non-Precious

Sales Representative: